




PROCAN ELECTRONICS INC.



PE - 6800 VET AUTO HEMATOLOGY ANALYZER

- **Methodology:** Electrical resistance for counting, hemiglobincyanide method and SFT method for hemoglobin
- **Parameter:** 3-part differentiation of WBC, 20 parameters and 3 color histograms(WBC,RBC,PLT)
- **Sample volume:** 9.6 L for venous mode, 20 L for prediluted mode
- **Species of animals:** Dogs, Cats, Horses, Pigs, Cattles, Goats, Monkey, Mouse etc. Additionally, four user-defined test modes are provided to meet the needs of various users
- **Quality Control:** L-j curve, X SD CV

[Http://www.procan.com.cn](http://www.procan.com.cn)

PROCAN